

Comparative Post Translational Modification in mAbdi2 and Innovator samples by LC-MS/MS workflow



Shrikant Mane^a, Şefik Önder^b, Halil İbrahim Özdemir^{a,c}, Aykut Demirkiran^a, Ravi Kumar Lella^a

^aAbdi İbrahim İlaç San. A.Ş. Orhan Gazi District, Tunç Street, No: 3 Esenyurt – İstanbul

^bYıldız Technical University, İstanbul

^cMarmara University, Bioengineering Department, Maltepe – İstanbul

shrikant.mane@abdiibrahim.com.tr

Introduction

The mAbdi2 is being developed as a biosimilar candidate. This study was designed to provide a comparative analytical assessment of the mAbdi2 and Innovator products to evaluate the post-translational modifications by LC-MS/MS workflow.

Methodology

The Bruker's maXis II-QTOF mass spectrometer was used in combination with an Elute UPLC in the positive ion polarity mode. Spectra were recorded in the mass range between 500-5000 m/z for Intact samples, and 300-3000 m/z for reduced samples. The de-glycosylation of the samples was performed by PNGase-F for overnight incubation at 37 °C. The reduction of the samples was performed by 1M DTT, at 65 °C for 15 Minutes. The spectra were deconvoluted using the Maximum Entropy algorithm. The mobile phase A 0.1% (v/v) formic acid in water and mobile phase B 0.1% (v/v) formic acid in acetonitrile.

For peptide mapping, samples were denatured with guanidine hydrochloride and then treated with dithiothreitol (DTT) to break disulfide linkages between cysteine residues for complete reduction. To prevent the disulfide bonds from re-forming, iodoacetamide (IAA) was used to modify the reactive cysteine-SH groups, forming carbamidomethylated cysteines. After alkylation, buffer exchange was performed for the removal of excess salts and trypsin was used for digesting the protein into peptides and identified by the mass spectrometer by Biopharma compass software.









Figure 1. LC-MS/MS workflow

Results



Post Translational Modifications (PTMS, site and extent)

- The PMF of the samples under Glyco & De-glyco reduced conditions were performed for the PTM analysis.
- The Oxidation, Deamidation, and Glycosylation site and extent of modifications were evaluated.
- The comparative N and C terminal modifications in the Light & Heavy chain for all the tested samples were presented.



- Similar PTMs observed in mAbdi2 & Innovator samples.

N and C Terminal Modifications (Heavy Chain)

- The mAbdi2 samples show comparatively less percentage of C-Term Lysine truncation.
- Whereas the relatively higher conversion of Glutamine (Q) to Pyro Glutamine (Pyro Q) in mAbdi2 was observed.
- For further confirmation, all the samples were treated with Carboxypeptidase B (CpB) & Glutaminyl-peptide cyclotransferase (QPCT) enzymes
- After treatment, samples were evaluated for charge variant profile and summarized below




Figure 5. Charge variants chromatograms with & without QPCT and Cp-B treatment

N-Term Q/E to Pyro Q Conversion v/s Relative Potency

- To understand the conversion mechanism of Q/E to Pyro Q under physiological conditions, a PBS incubation study was conducted.
- The samples were buffer exchanged to PBS, pH 7.4, and incubated at 37°C for 2 weeks.
- The CEX-HPLC profile was evaluated at 0, 1, and 2-week time points.
- The same PBS incubation study samples were used to perform the relative potency by Reporter Gene Assay and results were summarized.




Figure 7. Comparative Relative Potency by Reporter Gene Assay

Figure 6. CEX-HPLC profile under physiological condition (PBS-Incubation study)

Figure 7. Comparative Relative Potency by Reporter Gene Assay

PBS Incubation Time

Relative Potency % (mean n=3)

Samples

mAbdi2 Innovator US

Acceptance Criteria 70-130%

Figure 7. Comparative Relative Potency by Reporter Gene Assay

PBS Incubation Time

Relative Potency % (mean n=3)

Samples

mAbdi2 Innovator US

Acceptance Criteria 70-130%